

h5pydantic, from a Synchrotron through Python to HDF5

Clinton Roy

ANSTO
Australian Synchrotron

PyCon AU 2023

Synchrotron


HDF5

h5pydantic

Australian Synchrotron



Accelerating and bending electrons



- 1 Electron Gun
- 2 Linear Accelerator
- 3 Booster Ring
- 4 Storage Ring
- 5 Beamline
- 6 End Station

Lots of [meta]data

5

- ▶ sample information
 - ▶ chemistry, volume, concentration, thickness
 - ▶ strain, temperature, shearing, magnetic, UV light, chemical, pressure

Lots of [meta]data

- ▶ sample information
 - ▶ chemistry, volume, concentration, thickness
 - ▶ strain, temperature, shearing, magnetic, UV light, chemical, pressure
- ▶ x-ray beam setup
 - ▶ energy
 - ▶ distance between sample and detector
 - ▶ xray detector images
 - ▶ detector exposure time

Lots of [meta]data

- ▶ sample information
 - ▶ chemistry, volume, concentration, thickness
 - ▶ strain, temperature, shearing, magnetic, UV light, chemical, pressure
- ▶ x-ray beam setup
 - ▶ energy
 - ▶ distance between sample and detector
 - ▶ xray detector images
 - ▶ detector exposure time
- ▶ calibration
 - ▶ detector images with no beam
 - ▶ detector images with beam, but no sample
 - ▶ position of beam on detector

- ▶ sample information
 - ▶ chemistry, volume, concentration, thickness
 - ▶ strain, temperature, shearing, magnetic, UV light, chemical, pressure
- ▶ x-ray beam setup
 - ▶ energy
 - ▶ distance between sample and detector
 - ▶ xray detector images
 - ▶ detector exposure time
- ▶ calibration
 - ▶ detector images with no beam
 - ▶ detector images with beam, but no sample
 - ▶ position of beam on detector
- ▶ processing arguments
 - ▶ different statistical methods and arguments
 - ▶ different assumptions

- ▶ All the data and metadata together
- ▶ Ideally in one file

- ▶ One file, with internal structure
- ▶ Archive file, think tar/zip

1. Groups

- ▶ which are basically directories/folders

1. Groups
 - ▶ which are basically directories/folders
2. Datasets
 - ▶ multi dimensional arrays
 - ▶ 2D arrays commonly used to store images

1. Groups
 - ▶ which are basically directories/folders
2. Datasets
 - ▶ multi dimensional arrays
 - ▶ 2D arrays commonly used to store images
3. Attributes
 - ▶ key, value metadata
 - ▶ can attach to Groups and Datasets

```
/sample
/sample/thickness = 2.5
/sample/concentration = 3.0
/beam
/beam/energy = 1.8
/beam/distance = 7.8
/calibration
/calibration/no_beam = [0 0 0 0 1 0 0 1 0 ...]
/calibration/position/x = 247
/calibration/position/y = 253
/processing
/processing/threshold = 25
```

```
import h5py

h5file = h5py.File("experiment.hdf5", "w")

calibration_group = h5file.create_group("calibration")

dataset = calibration_group.create_dataset("no_beam",
                                           [0, 0, 0, 1, 0, 0])

position = calibration_group.create_group("position")
position.attrs["x"] = 247
position.attrs["y"] = 253
```


- ▶ Groups - Model
- ▶ Dataset - Model
- ▶ Attributes - fields of Groups/Datasets

```
from h5pydantic import H5Group, H5Dataset, H5Integer32

class Position(H5Group):
    x: H5Integer32
    y: H5Integer32

class CalibImage(H5Dataset, dtype=H5Integer32, shape=(3,2)):
    pass

class Calibration(H5Group):
    no_beam: CalibImage = CalibImage()
    position: Position

class Experiment(H5Group):
    calibration: Calibration
```

```
from model import Experiment

exp = Experiment(calibration={"position": {"x": 247, "y": 253}})
exp.calibration.no_beam.data([[0, 0, 1], [0, 0, 0]])
exp.dump("experiment.hdf5")
```

And later on...

```
from model import Experiment

exp = Experiment.load("experiment.hdf5")
(x, y) = exp.calibration.position.x, exp.calibration.position.y
```

- ▶ Can extend pydantic fields with all sorts of information
 - ▶ Documentation
 - ▶ Units
 - ▶ Relationships between fields

```
from h5pydantic import H5Group, H5Integer32
from pydantic import Field

import pint
ureg = pint.UnitRegistry()

class Beam(H5Group):
    energy: float = Field(ge=0, doc="X-Ray beam energy", unit=ureg.joule)
    distance: H5Integer32 = Field(ge=0,
                                  doc="Distance from sample to detector",
                                  unit=ureg.millimeter)

beam = Beam(energy=3.4, distance=7500)
```

- ▶ Can output a HDF5 description in PDF, using Sphinx
- ▶ Can start to handle versioning

The overall design of `h5pydantic` has three main aims:

1. All instantiated models can be saved to HDF5 files
2. All HDF5 features can be used.
3. All Pydantic features can be used.

Questions?

- ▶ Thank you to the scientists that have test driven the library:
Dr Lester Barnsley
- ▶ Thank you to all the members of the Scientific Computing team and other Scientists for feedback on the library and this presentation